




# Color Signal and Spectrometry

## Light Source

Spectral Exitance

$M_{e,\lambda}$  in  $\text{W}\cdot\text{m}^{-2}\cdot\text{nm}^{-1}$   
("Emittance")




## Object

Spectral Reflectance  $R_\lambda$  in % or (0 to 1)

## Color Signal or Spectral Radiosity

$$J_{e,\lambda} = M_{e,\lambda} R_\lambda$$

Also in  $\text{W}\cdot\text{m}^{-2}\cdot\text{nm}^{-1}$




Thus, to get Reflectance


$$R_{\lambda} = \frac{J_{e,\lambda}}{M_{e,\lambda}}$$


Spectrometers can automatically do this calculation.  
But first you have to measure emittance and dark current.

Pointing the spectrometer straight to the light source can saturate its sensor. Instead a white reference is used.




# Lambertian Surface - Ideal matte surface





# To get reflectance using a white reference standard



NOTE: Spectral Radiance is a DIRECTIONAL quantity.


To get reflectance of object

$$R_{\lambda} = \frac{L(object)_{e,\Omega,\lambda}}{L(white)_{e,\Omega,\lambda}}$$




# Use of Integrating Sphere for Spectrometry

## 1. Measuring transmittance



## 2. Measuring reflectance



- Instead of **integrating sphere** you will use your **lightbox**.
- Instead of a **spectralon reference** you will use several layers of **white matte paper**.
- To remove directional dependence, **fix the optical fiber on a holder**.